

SMART GLOVE™

**Product
Catalogue**



Your Futuristic Gloves



Smart Glove Corporation



Sigma Glove and GX3 Specialty Plant



GX Corporation and Platinum Glove

About SMART GLOVE™

ABOUT US

Founded in 1994, Smart Glove is one of Malaysia's leading innovative glove manufacturers recognised on a global front. Spearheaded by the founder of Smart Glove, Mr KP Foo leads his team of experts with his sound technical knowledge and hands-on approach in Research & Development. Armed with a wealth of experience and expertise in glove manufacturing, and with a robust operating history of over 20 years, Smart Glove provides a diverse range of gloves with different specifications to meet the needs of our customers from various industries including medical, food and beverage, and cleanroom among others.

Smart Glove Group is currently operating in Malaysia and Indonesia, which are operationalized with modern purpose-built facilities.

OUR MISSION

As the pioneer of Nitrile glove technology, Smart Glove strives to be at the forefront of the glove industry by continuously researching and developing gloves that are reliable, comfortable, environmentally friendly and suitable for various applications.

OUR COMMITMENTS

Our gloves undergo stringent quality inspections and independent third party assessments prior to each shipment. This indicates our commitment to continuously exceed the Acceptable Quality Level (AQL) standard, providing our customers with maximum barrier protection and consistently delivering products that surpasses customers' quality requirements.



OUR PARTNERS

Our partners have trusted us as their key OEM manufacturer and single-use hand protection solutions provider for over two decades.

All the companies in the group are ISO certified, embracing good manufacturing practices, be it with products meeting or surpassing various international standards, such as The Food and Drug Administration (FDA) and European Directives (EN), as well as meeting market regulations in Japan, Canada, Brazil, Russia, and many more.



Smart Glove is committed to ensuring and delivering the highest quality of its gloves and has been accredited with various internationally recognised certifications.

International Certification

Food Contact Regulation

EN ISO 374-1:2016/type C





Medical

Our exam gloves are equipped with a new improved technology, designed to give superior comfort and easy donning

- **Suitable for** medical personnel, veterinary operations, pharmaceutical and general hygiene use
- **Offers a variety of colours** to choose from
- **Investigated and verified by UL** in accordance with the Standard on Protective Clothing and Ensembles for Emergency Medical Operations, NFPA 1999-2018

Material

Nitrile, Polychloroprene

Length

9", 12"

Colour

Black, White, Pink, Green, Dark Green, Blue, Steel Blue, Teal Blue, Dark Blue, Cobalt Blue, Ice Blue, Violet Blue

Thickness

3mil to 10mil

Dark Green 9"
Polychloroprene C-Series
glove



Blue 9"
Nitrile S-Series glove





Dental

The classic dental gloves are safe for oral examination to deliver a valuable patient experience

- **Flexible gloves** give ultimate wet grip and comfort for best performance
- **Available in a variety of scents** to choose from
- **100% Nitrile glove.** Safe for use from latex protein allergy

Material

Nitrile, Polychloroprene

Length

9"

Colour

Blue, Violet Blue, Pink, White, Orange, Green

Scent

Grape, Mint, Bubble Gum, Green tea, Vanilla, Orange

Thickness

3mil to 6mil






Surgical

Our surgical gloves are designed and custom-made to meet the high expectations of surgeons

- **New improved technology** designed to give superior softness to aid easy donning and to reduce hand fatigue
- **We offer various types of materials** in order to meet the different requirements for different types of surgeries
- **Sterile Exam gloves** are also available and can be packed in singles or pairs per pouch
- **Chemotherapy Drugs Resistance** - our gloves are tested according to the Chemotherapy Permeation test



Surgical

Material
Nitrile, Polychloroprene, Polyisoprene, Natural Rubber

Length
12"

Colour
White, Cobalt Blue, Dark Blue, Green, Dark Green, Brown

Thickness
8mil, 6mil (Underglove)

Sterile Exam

Material
Nitrile, Polychloroprene

Length
9", 12"

Colour
Blue, White

Thickness
3.5mil, 6mil





Laboratory / EMS

Our outstanding gloves are designed with chemical resistance for laboratory and Emergency Medical Services use

- **Gloves dipped with dual colour** provide enhanced protection for better visibility of any cuts or pricks when in contact with fluid or blood
- **Provides superior strength** to ensure optimal barrier protection
- **Offers hybrid material dipped** to provide chemical permeation resistance

Material

Nitrile, Polychloroprene, Natural Rubber, Hybrid material (Polychloroprene (exterior) Nitrile (interior))

Length

9", 12", 16"

Colour

Violet Blue, Blue
Dual Colour: Violet Blue / Green, Black / Red, Orange / Black, Black / Silver

Thickness

3mil to 18mil
5mil to 13mil (Hybrid material)
15mil to 18mil (Natural Rubber)





High Risk

The ideal glove designed for high risk exposure and heavy duty tasks

- **Featuring Fish Scale and Diamond textured gloves** specially designed with anti-slip for excellent traction during wet and oily applications
- **Suitable for** automotive repair and maintenance, spray paint, machinery maintenance and chemical handling
- **Offers extra glove thickness of up to 20mil** for better protection

Material

Nitrile, Polychloroprene, Natural Rubber

Length

9", 12", 16"

Colour

Black, Orange, Dark Blue, Green, Dual Colour


Thickness

6mil to 20mil

Dark Blue Fish Scale 9"
Nitrile glove



Orange Diamond 12"
Nitrile glove





Cleanroom

Our cleanroom gloves feature exceptionally low levels of particles as well as ionic extractable and non-volatile residue

Material
Nitrile

Length
240mm, 300mm, 330mm

Colour
Opaque White, Translucent

Thickness
3.5mil, 4.0mil, 4.5mil, 6-7mil

- **New formulation** designed to provide comfort when executing tasks over long hours
- **Suitable for** Semiconductor, Microelectronic, Aerospace, Optics and Biotechnology, Hard Disk Drive, Healthcare, Laboratory, Pharmaceutical and manufacturing industries that require a clean and controlled environment
- **Options of** Class 10, Class 100 and Class 1000
- Offers superior **electrostatic dissipative (ESD) properties** (Antistatic)





Food

Our gloves are suitable for handling all types of food

- **Food Processing** gloves are designed to provide strong grip in handling wet, dry and oily foods. Suitable for fish, chicken and poultry processing
- Available in **Fish Scale** and **Diamond** textures
- **Food Service** gloves are designed to provide excellent dexterity and performance. Suitable for cooking, cutting, slicing and preparing food

Food Processing

Material
Nitrile, Polychloroprene, Natural Rubber

Length
9", 12", 16"

Colour
Black, Orange, Dark Blue, Dual Colour, Dark Steel Blue

Thickness
6mil to 10mil

Food Service

Material
Nitrile

Length
9"

Colour
White, Dark Steel Blue, Translucent, Smart Blue

Thickness
2mil

Black Fish Scale 9"
Nitrile glove



Translucent & Dark Steel Blue 9"
Nitrile glove



Pattern

The choices are yours - from colour, material, size to texture.
We can customise according to your requirements.

For more information, please contact us at info@smartglove.com.my

Finger or Full Textured

Fish Scale

Diamond


Mix and match to make it your own

Choosing the right gloves for you.



1. Material

- Nitrile
- Polychloroprene
- Polyisoprene
- Natural Rubber
- Hybrid
- Special Formulation



2. Colour

- Various Solid Colours
- Dual Colours



3. Pattern

- Finger Textured
- Full Textured
- Fish Scale Textured
- Diamond Textured
- Micro-roughened



4. Others

- Sterile
- Extended Cuff
- Scented
- Hand Specific
- Chlorinated or Polymer Coated



SMART GLOVE CORPORATION SDN BHD

Lot 6487, Batu 5 ¼ Sementa, Jalan Kapar,
42100 Klang, Selangor, Malaysia

Tel : +603 - 3291 6100 Email : info@smartglove.com.my
Fax : +603 - 3291 7526 www.smartglove.com.my

*Smart Glove is a trademark owned by Smart Glove Holdings Sdn Bhd or one of its affiliates.
© 2018 Smart Glove Holdings Sdn Bhd. All Rights Reserved.*