




# SmartMedic™



*Our surgical gloves are designed and custom-made to meet the high expectations of surgeons. We offer various types of materials in order to meet the different requirement for different types of surgeries.*

»»»» Provides superior softness to aid easy donning and to reduce hand fatigue

»»»» Sterile Exam gloves are available in singles or pairs per pouch packaging

»»»» Tested safe for use with chemotherapy drugs per global Chemotherapy Permeation Test standard



**SmartMedic™**  
Sterile Surgical Glove



<b>Material</b>	Nitrile/ Polychloroprene/ Polyisoprene/ Natural Rubber
<b>Length</b>	12-inch
<b>Finishing</b>	Powdered/ Powder Free
<b>Design &amp; Features</b>	Hand Specific, Micro Textured, Beaded Cuff
<b>Colours</b>	Blue/ White/ Green
<b>Quality Standards</b>	ASTM D3577, ASTM D6978, EN455 1-4, EN420
<b>Weight (M Size)</b>	8.6g - 12.5g (6mil underglove - 8mil)
<b>Highlights</b>	The hand specific gloves are designed considering the anatomy of human hand in order to provide more comfort and less fatigue under prolonged time of wearing Damp Hand Don feature available Among the various surgical gloves, the surgical glove made out of Polyisoprene is the most comfortable to the surgeon

**SmartMedic™**  
Sterile Examination Glove



<b>Material</b>	Nitrile
<b>Length</b>	9-inch, 12-inch
<b>Finishing</b>	Powder Free
<b>Design &amp; Features</b>	Ambidextrous/ Finger Textured or Full Textured / Beaded Cuff
<b>Colours</b>	Blue/ White
<b>Quality Standards</b>	ASTM D6319, ASTM D6978, EN455 1-4, EN420
<b>Weight (M Size)</b>	4.0g - 6.6g (3.5mil - 6mil)



**Standards & Compliances**

**ASTM D3577**

Standard Specification for Rubber Surgical Gloves

**ASTM D6978-05 (R2019)**

Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs

**EN 455**

Requirements and Test for Single use medical gloves

**Regulation (EU) 2016/425**

Personal Protective Equipment

**EN ISO 21420**

Protective Gloves: General Requirements and Test Methods

\*Other specific local regulatory requirements can be fulfilled upon request.

**Why SmartMedic™**



Excellent sensitivity and grip



Superior comfort and easy donning to the wearer



Advanced allergy protection



Protection against contamination

**Applications**



Surgical Procedure



Cleanroom environment



Controlled environment

Our Brands:

**GEN-X®**   **iDentalglove™**   **SmartMedic™**

**MetalMan™**   **iChefglove™**   **SmartClean™**   **iSafe™**